

## Sialin

Catalog No: #AP87592



Package Size: #AP87592-1 50ug #AP87592-2 100ug #AP87592-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sialin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	H(+)/nitrate cotransporter,H(+)/sialic acid cotransporter,AST,Membrane glycoprotein HP59,Solute carrier family 17 member 5,Vesicular H(+)/Aspartate-glutamate cotransporter
Accession No.	Uniprot ID:Q9NRA2
Uniprot	Q9NRA2
GeneID	26503;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SLC17A5

Note: This product is for in vitro research use only