

## Glutaminyl-peptide cyclotransferase-like protein

Catalog No: #AP87599



Package Size: #AP87599-1 50ug #AP87599-2 100ug #AP87599-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glutaminyl-peptide cyclotransferase-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Golgi-resident glutaminyl-peptide cyclotransferase,isoQC,gQC
Accession No.	Uniprot ID:Q0V8G3
Uniprot	Q0V8G3
GeneID	530871;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:QPCTL

Note: This product is for in vitro research use only