

Ras-related protein Rap-1A

Catalog No: #AP87601



Package Size: #AP87601-1 50ug #AP87601-2 100ug #AP87601-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Ras-related protein Rap-1A |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | C21KG,G-22K,GTP-binding protein smg p21A,Ras-related protein Krev-1,KREV1 |
| Accession No. | Uniprot ID:P62834 |
| Uniprot | P62834 |
| GeneID | 5906; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RAP1A

Note: This product is for in vitro research use only