

Retinoic acid receptor gamma

Catalog No: #AP87602



Package Size: #AP87602-1 50ug #AP87602-2 100ug #AP87602-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	RAR-gamma,Nuclear receptor subfamily 1 group B member 3,NR1B3
Accession No.	Uniprot ID:P13631
Uniprot	P13631
GeneID	5916;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RARG

Note: This product is for in vitro research use only