

## Serpín H1

Catalog No: #AP87608

Package Size: #AP87608-1 50ug #AP87608-2 100ug #AP87608-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serpín H1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	47 kDa heat shock protein, Collagen-binding protein, Colligin, HSP47
Accession No.	Uniprot ID:P13731
Uniprot	P13731
GeneID	396228;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: SERPINH1

Note: This product is for in vitro research use only