

Sigma non-opioid intracellular receptor 1

Catalog No: #AP87609



Package Size: #AP87609-1 50ug #AP87609-2 100ug #AP87609-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sigma non-opioid intracellular receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Sigma 1-type opioid receptor, Sigma1-receptor, OPRS1, RCJMB04_7e2
Accession No.	Uniprot ID:Q5ZL84
Uniprot	Q5ZL84
GeneID	100859748;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SIGMAR1

Note: This product is for in vitro research use only