

Regucalcin

Catalog No: #AP87610



Package Size: #AP87610-1 50ug #AP87610-2 100ug #AP87610-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Regucalcin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	RC,Gluconolactonase,GNL,Senescence marker protein 30,SMP-30,SMP30
Accession No.	Uniprot ID:Q9TTJ5
Uniprot	Q9TTJ5
GeneID	280910;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RGN

Note: This product is for in vitro research use only