

Sperm-associated antigen 1

Catalog No: #AP87616



Package Size: #AP87616-1 50ug #AP87616-2 100ug #AP87616-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sperm-associated antigen 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Infertility-related sperm protein Spag-1,TPR-containing protein involved in spermatogenesis,TPIS,Tpis
Accession No.	Uniprot ID:Q80ZX8
Uniprot	Q80ZX8
GeneID	26942;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Spag1

Note: This product is for in vitro research use only