

Mothers against decapentaplegic homolog 5

Catalog No: #AP87617



Package Size: #AP87617-1 50ug #AP87617-2 100ug #AP87617-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MAD homolog 5, SMAD family member 5, SMAD 5, Madh5
Accession No.	Uniprot ID:Q9R1V3
Uniprot	Q9R1V3
GeneID	59328;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Smad5

Note: This product is for in vitro research use only