

Reactive oxygen species modulator 1

Catalog No: #AP87623



Package Size: #AP87623-1 50ug #AP87623-2 100ug #AP87623-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Reactive oxygen species modulator 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	ROS modulator 1,Protein MGR2 homolog
Accession No.	Uniprot ID:A1XQR6
Uniprot	A1XQR6
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ROMO1

Note: This product is for in vitro research use only