

Inactive tyrosine-protein kinase transmembrane receptor ROR1

Catalog No: #AP87624

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87624-1 50ug #AP87624-2 100ug #AP87624-3 1mg

Description

Product Name	Inactive tyrosine-protein kinase transmembrane receptor ROR1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	mROR1, Neurotrophic tyrosine kinase, receptor-related 1, Ntrkr1
Accession No.	Uniprot ID: Q9Z139
Uniprot	Q9Z139
GeneID	26563;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ror1

Note: This product is for in vitro research use only