

Protein S100-A2

Catalog No: #AP87625



Package Size: #AP87625-1 50ug #AP87625-2 100ug #AP87625-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein S100-A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Protein S-100L,S100 calcium-binding protein A2
Accession No.	Uniprot ID:P10462
Uniprot	P10462
GeneID	509860;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S100A2

Note: This product is for in vitro research use only