

Protein S100-A6

Catalog No: #AP87626



Package Size: #AP87626-1 50ug #AP87626-2 100ug #AP87626-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein S100-A6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Calcyclin, Prolactin receptor-associated protein, S100 calcium-binding protein A6, Cacy
Accession No.	Uniprot ID:P05964
Uniprot	P05964
GeneID	85247;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: S100a6

Note: This product is for in vitro research use only