

## Bile acyl-CoA synthetase

Catalog No: #AP87628



Package Size: #AP87628-1 50ug #AP87628-2 100ug #AP87628-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bile acyl-CoA synthetase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BACS,Bile acid-CoA ligase,BA-CoA ligase,Cholate--CoA ligase,Fatty acid transport protein 5,FATP-5,Fatty-acid-coenzyme A ligase, very long-chain 3,Solute carrier family 27 member 5,Very long-chain acyl-CoA synthetase homolog 2,VLCS-H2,Very long-chain acyl-CoA synthetase-related protein,VLACS-related,ACSB,ACSVL6,FACVL3,FATP5
Accession No.	Uniprot ID:Q9Y2P5
Uniprot	Q9Y2P5
GeneID	10998;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SLC27A5

Note: This product is for in vitro research use only