

## Sodium channel protein type 9 subunit alpha

Catalog No: #AP87631



Package Size: #AP87631-1 50ug #AP87631-2 100ug #AP87631-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sodium channel protein type 9 subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Peripheral sodium channel 1,PN1,Sodium channel protein type IX subunit alpha,Voltage-gated sodium channel subunit alpha Nav1.7,Kiaa4197
Accession No.	Uniprot ID:Q62205
Uniprot	Q62205
GeneID	20274;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Scn9a

Note: This product is for in vitro research use only