

Radical S-adenosyl methionine domain-containing protein 2

Catalog No: #AP87632

Package Size: #AP87632-1 50ug #AP87632-2 100ug #AP87632-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Radical S-adenosyl methionine domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Bone-expressed sequence tag 5 protein,Viperin,Virus inhibitory protein, endoplasmic reticulum-associated, interferon-inducible,Best5
Accession No.	Uniprot ID:O70600
Uniprot	O70600
GeneID	65190;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rsad2

Note: This product is for in vitro research use only