

## Retinoic acid receptor RXR-beta

Catalog No: #AP87633



Package Size: #AP87633-1 50ug #AP87633-2 100ug #AP87633-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Retinoic acid receptor RXR-beta  |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Hu   |
| Other Names        | Nuclear receptor subfamily 2 group B member 2, Retinoid X receptor beta, NR2B2 |
| Accession No.      | Uniprot ID:P28702  |
| Uniprot            | P28702   |
| GeneID             | 6257;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                         |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: RXRB

Note: This product is for in vitro research use only