

Proton-associated sugar transporter A

Catalog No: #AP87634



Package Size: #AP87634-1 50ug #AP87634-2 100ug #AP87634-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proton-associated sugar transporter A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PAST-A, Deleted in neuroblastoma 5 protein, DNB-5, Solute carrier family 45 member 1, DNB5
Accession No.	Uniprot ID:Q9Y2W3
Uniprot	Q9Y2W3
GeneID	50651;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC45A1

Note: This product is for in vitro research use only