

Mothers against decapentaplegic homolog 3

Catalog No: #AP87635



Package Size: #AP87635-1 50ug #AP87635-2 100ug #AP87635-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	MAD homolog 3, SMAD family member 3, SMAD 3, MADH3
Accession No.	Uniprot ID:P84024
Uniprot	P84024
GeneID	397260;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SMAD3

Note: This product is for in vitro research use only