

Serotonin N-acetyltransferase

Catalog No: #AP87636



Package Size: #AP87636-1 50ug #AP87636-2 100ug #AP87636-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serotonin N-acetyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Serotonin acetylase,Aralkylamine N-acetyltransferase,AA-NAT,Snat
Accession No.	Uniprot ID:Q64666
Uniprot	Q64666
GeneID	25120;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Aanat

Note: This product is for in vitro research use only