

Signal-induced proliferation-associated 1-like protein 1

Catalog No: #AP87637



Package Size: #AP87637-1 50ug #AP87637-2 100ug #AP87637-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Signal-induced proliferation-associated 1-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SIPA1-like protein 1,Kiaa0440
Accession No.	Uniprot ID:Q8C0T5
Uniprot	Q8C0T5
GeneID	217692;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sipa111

Note: This product is for in vitro research use only