

Semaphorin-3E

Catalog No: #AP87640



Package Size: #AP87640-1 50ug #AP87640-2 100ug #AP87640-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Semaphorin-3E
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Collapsin-5, COLL-5, COLL5
Accession No.	Uniprot ID:O42237
Uniprot	O42237
GeneID	374089;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SEMA3E

Note: This product is for in vitro research use only