

Semaphorin-3F

Catalog No: #AP87641



Package Size: #AP87641-1 50ug #AP87641-2 100ug #AP87641-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Semaphorin-3F
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Semaphorin IV, Sema IV
Accession No.	Uniprot ID:O88632
Uniprot	O88632
GeneID	20350;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sema3f

Note: This product is for in vitro research use only