

Septin-5

Catalog No: #AP87642



Package Size: #AP87642-1 50ug #AP87642-2 100ug #AP87642-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Septin-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cell division control-related protein 1,CDCrel-1,Peanut-like protein 1,Pnutl1
Accession No.	Uniprot ID:Q9Z2Q6
Uniprot	Q9Z2Q6
GeneID	18951;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:43713

Note: This product is for in vitro research use only