

Secretory carrier-associated membrane protein 1

Catalog No: #AP87644



Package Size: #AP87644-1 50ug #AP87644-2 100ug #AP87644-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Secretory carrier-associated membrane protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Secretory carrier membrane protein 1
Accession No.	Uniprot ID:O77735
Uniprot	O77735
GeneID	397477;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SCAMP1

Note: This product is for in vitro research use only