

Ribosome-binding protein 1

Catalog No: #AP87647



Package Size: #AP87647-1 50ug #AP87647-2 100ug #AP87647-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ribosome-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	180 kDa ribosome receptor,RRp,P180
Accession No.	Uniprot ID:Q28298
Uniprot	Q28298
GeneID	403809;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RRBP1

Note: This product is for in vitro research use only