

Adenosylhomocysteinase

Catalog No: #AP87650



Package Size: #AP87650-1 50ug #AP87650-2 100ug #AP87650-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Adenosylhomocysteinase |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | AdoHcyase,S-adenosyl-L-homocysteine hydrolase |
| Accession No. | Uniprot ID:Q3MHL4 |
| Uniprot | Q3MHL4 |
| GeneID | 508158; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AHCY

Note: This product is for in vitro research use only