

# Endophilin-A2

Catalog No: #AP87652



Package Size: #AP87652-1 50ug #AP87652-2 100ug #AP87652-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Endophilin-A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Endophilin-2,SH3 domain protein 2B,SH3 domain-containing GRB2-like protein 1,SH3p8,Sh3p8
Accession No.	Uniprot ID:O35964
Uniprot	O35964
GeneID	81922;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Sh3gl1

Note: This product is for in vitro research use only