

Dexamethasone-induced Ras-related protein 1

Catalog No: #AP87654



Package Size: #AP87654-1 50ug #AP87654-2 100ug #AP87654-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dexamethasone-induced Ras-related protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Activator of G-protein signaling 1,AGS1,DEXRAS1
Accession No.	Uniprot ID:Q9Y272
Uniprot	Q9Y272
GeneID	51655;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RASD1

Note: This product is for in vitro research use only