

Retinoblastoma-like protein 2

Catalog No: #AP87656



Package Size: #AP87656-1 50ug #AP87656-2 100ug #AP87656-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinoblastoma-like protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	130 kDa retinoblastoma-associated protein,p130,Retinoblastoma-related protein 2,RBR-2,pRb2,RB2
Accession No.	Uniprot ID:Q08999
Uniprot	Q08999
GeneID	5934;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RBL2

Note: This product is for in vitro research use only