

# Semaphorin-4D

Catalog No: #AP87658



Package Size: #AP87658-1 50ug #AP87658-2 100ug #AP87658-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Semaphorin-4D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	M-Sema G,Semaphorin-C-like 2,Semaphorin-J,Sema J,CD100,Semacl2,Semaj
Accession No.	Uniprot ID:O09126
Uniprot	O09126
GeneID	20354;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Sema4d

Note: This product is for in vitro research use only