

## Double-stranded RNA-specific editase 1

Catalog No: #AP87661



Package Size: #AP87661-1 50ug #AP87661-2 100ug #AP87661-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Double-stranded RNA-specific editase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	RNA-editing deaminase 1, RNA-editing enzyme 1, dsRNA adenosine deaminase, Red1
Accession No.	Uniprot ID:P51400
Uniprot	P51400
GeneID	25367;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Adarb1

Note: This product is for in vitro research use only