

# Protein S100-A11

Catalog No: #AP87664



Package Size: #AP87664-1 50ug #AP87664-2 100ug #AP87664-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein S100-A11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Calgizzarin, Protein S100-C, S100 calcium-binding protein A11, S100C
Accession No.	Uniprot ID:P31950
Uniprot	P31950
GeneID	445534;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: S100A11

Note: This product is for in vitro research use only