

Protein S100-A12

Catalog No: #AP87665



Package Size: #AP87665-1 50ug #AP87665-2 100ug #AP87665-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein S100-A12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CGRP, Calcium-binding protein in amniotic fluid 1, CAAF1, Calgranulin-C, CAGC, Extracellular newly identified RAGE-binding protein, EN-RAGE, Migration inhibitory factor-related protein 6, MRP-6, Neutrophil S100 protein, S100 calcium-binding protein A12
Accession No.	Uniprot ID:P80511
Uniprot	P80511
GeneID	6283;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: S100A12

Note: This product is for in vitro research use only