

Secretogranin-3

Catalog No: #AP87668



Package Size: #AP87668-1 50ug #AP87668-2 100ug #AP87668-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Secretogranin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Secretogranin III, SgIII
Accession No.	Uniprot ID:A6QLI2
Uniprot	A6QLI2
GeneID	525941;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SCG3

Note: This product is for in vitro research use only