

NAD-dependent deacetylase sirtuin-6

Catalog No: #AP87669



Package Size: #AP87669-1 50ug #AP87669-2 100ug #AP87669-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD-dependent deacetylase sirtuin-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SIR2-like protein 6,SIR2L6
Accession No.	Uniprot ID:Q8N6T7
Uniprot	Q8N6T7
GeneID	51548;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SIRT6

Note: This product is for in vitro research use only