

Retinoic acid receptor RXR-alpha

Catalog No: #AP87671



Package Size: #AP87671-1 50ug #AP87671-2 100ug #AP87671-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor RXR-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Nuclear receptor subfamily 2 group B member 1, Retinoid X receptor alpha, Nr2b1
Accession No.	Uniprot ID:Q05343
Uniprot	Q05343
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Rxra

Note: This product is for in vitro research use only