

Protein S100-A10

Catalog No: #AP87673



Package Size: #AP87673-1 50ug #AP87673-2 100ug #AP87673-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein S100-A10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Calpactin I light chain, Calpactin-1 light chain, Cellular ligand of annexin II, Nerve growth factor-induced protein 42C, S100 calcium-binding protein A10, p10 protein, p11
Accession No.	Uniprot ID:P05943
Uniprot	P05943
GeneID	81778;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: S100a10

Note: This product is for in vitro research use only