

Retinoic acid receptor responder protein 2

Catalog No: #AP87678



Package Size: #AP87678-1 50ug #AP87678-2 100ug #AP87678-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor responder protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Chemerin,RAR-responsive protein TIG2,Tazarotene-induced gene 2 protein,TIG2
Accession No.	Uniprot ID:Q99969
Uniprot	Q99969
GeneID	5919;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RARRES2

Note: This product is for in vitro research use only