

## Serine/arginine-rich splicing factor 2

Catalog No: #AP87688



Package Size: #AP87688-1 50ug #AP87688-2 100ug #AP87688-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/arginine-rich splicing factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein PR264, Putative myelin regulatory factor 1, MRF-1, Splicing component, 35 kDa, Splicing factor SC35, SC-35, Splicing factor, arginine/serine-rich 2, Pr264, Sfrs10, Sfrs2
Accession No.	Uniprot ID:Q62093
Uniprot	Q62093
GeneID	20382;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Srsf2

Note: This product is for in vitro research use only