

Roundabout homolog 1

Catalog No: #AP87690

Package Size: #AP87690-1 50ug #AP87690-2 100ug #AP87690-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Roundabout homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Deleted in U twenty twenty,H-Robo-1,DUTT1
Accession No.	Uniprot ID:Q9Y6N7
Uniprot	Q9Y6N7
GeneID	6091;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ROBO1

Note: This product is for in vitro research use only