

Heterochromatin-associated protein MENT

Catalog No: #AP87695



Package Size: #AP87695-1 50ug #AP87695-2 100ug #AP87695-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heterochromatin-associated protein MENT
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Myeloid and erythroid nuclear termination stage-specific protein, Serpin B10, MENT
Accession No.	Uniprot ID:O73790
Uniprot	O73790
GeneID	395715;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SERPINB10

Note: This product is for in vitro research use only