

## Serpín B5

Catalog No: #AP87696



Package Size: #AP87696-1 50ug #AP87696-2 100ug #AP87696-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serpín B5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Maspín,Peptidase inhibitor 5,PI-5,PI5,Spi5,Spi7
Accession No.	Uniprot ID:P70124
Uniprot	P70124
GeneID	20724;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Serpínb5

Note: This product is for in vitro research use only