

Syndecan-1

Catalog No: #AP87699



Package Size: #AP87699-1 50ug #AP87699-2 100ug #AP87699-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Syndecan-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SYND1,CD138,Synd1
Accession No.	Uniprot ID:P26260
Uniprot	P26260
GeneID	25216;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sdc1

Note: This product is for in vitro research use only