

Protein S100-A7

Catalog No: #AP87701



Package Size: #AP87701-1 50ug #AP87701-2 100ug #AP87701-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein S100-A7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Calcium-binding protein in amniotic fluid 2,CAAF2,Dander minor allergen BDA11,Dermal allergen BDA11,S100 calcium-binding protein A7,CAAF2
Accession No.	Uniprot ID:Q28050
Uniprot	Q28050
GeneID	282344;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S100A7

Note: This product is for in vitro research use only