

Stathmin

Catalog No: #AP87702



Package Size: #AP87702-1 50ug #AP87702-2 100ug #AP87702-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Stathmin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Leukemia-associated gene protein,Leukemia-associated phosphoprotein p18,Metablastin,Oncoprotein 18,Op18,Phosphoprotein p19,pp19,Prosolin,Protein Pr22,pp17,Lag,Lap18,Pr22
Accession No.	Uniprot ID:P54227
Uniprot	P54227
GeneID	16765;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Stmn1

Note: This product is for in vitro research use only