

SH3 domain-containing YSC84-like protein 1

Catalog No: #AP87703



Package Size: #AP87703-1 50ug #AP87703-2 100ug #AP87703-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SH3 domain-containing YSC84-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:O08641
Uniprot	O08641
GeneID	24057;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sh3yl1

Note: This product is for in vitro research use only