

Serum amyloid A protein

Catalog No: #AP87704



Package Size: #AP87704-1 50ug #AP87704-2 100ug #AP87704-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serum amyloid A protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	SAA
Accession No.	Uniprot ID:P19708
Uniprot	P19708
GeneID	403585;751814;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SAA1

Note: This product is for in vitro research use only