

Mothers against decapentaplegic homolog 6

Catalog No: #AP87705



Package Size: #AP87705-1 50ug #AP87705-2 100ug #AP87705-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	MAD homolog 6, SMAD family member 6, SMAD 6, MADH6
Accession No.	Uniprot ID:Q9W734
Uniprot	Q9W734
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SMAD6

Note: This product is for in vitro research use only