

Serine/threonine-protein kinase SMG1

Catalog No: #AP87706



Package Size: #AP87706-1 50ug #AP87706-2 100ug #AP87706-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase SMG1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SMG-1,61E3.4,Lambda/iota protein kinase C-interacting protein,Lambda-interacting protein,ATX,KIAA0421,LIP
Accession No.	Uniprot ID:Q96Q15
Uniprot	Q96Q15
GeneID	23049;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SMG1

Note: This product is for in vitro research use only